

## Accessory for HAAKE MARS

### Universal Container Holder for HAAKE MARS and RheoStress 600

A universal holder for sample containers which can be used for the HAAKE MARS and the RheoStress 600 (as well as all of its predecessors) has been developed.

Thanks to the three individually adjustable clamps, sample containers like glass jars, cans, beakers, cups, etc., can be mounted solidly onto the rheometer. This allows for the insertion of a (vane) measuring geometry into the sample without comprising its structure.

Vane rotors with different diameters are available. The universal adapter shaft can be used to adapt any individually designed measuring geometry.

Using the HAAKE RheoWin rheometer software fully automatic measuring and evaluation routines (jobs), including quality check criteria can be created and run for fast quality control, e.g. yield point determination or viscosity measurement.

#### Order Information:

603-0304 Universal container holder for the HAAKE MARS and the HAAKE RheoStress 600 (as well as all of its predecessors). At least one measuring geometry is required.

#### Optional accessories:

222-1276 Adapter shaft with a 6 mm bore for vane, pin or individual rotors

222-1324 Vane rotor FL40 (D=40 mm)

222-1325 Vane rotor FL22 (D=22 mm)

222-1326 Vane rotor FL16 (D=16 mm)

222-1626 Rotor with variable number of laterally mounted pins

#### Key-words

- HAAKE MARS
- HAAKE RheoStress 600
- Universal Container

Dr. Cornelia Küchenmeister  
Jint Nijman

Thermo Fisher Scientific  
Process Instruments  
Dieselstr. 4  
76227 Karlsruhe  
Tel: +49 (0) 721 4 09 44 44  
Info.mc.de@thermofisher.com

[www.thermo.com/mc](http://www.thermo.com/mc)



Fig. 1: Universal container holder with three adjustable clamps (left),  
Fig. 2: Universal clamp holder with a peanut butter glass jar mounted on the HAAKE MARS with a vane rotor FL22 (right)