

Product:

GLACIER
Ultra Low Temperature Circulators

Manufactured By:

Thermo Fisher Scientific

| | G50-AC200 | G50-PC200 |
|---------------------------------|---------------------|---------------------|
| Temperature Range | -50°C - +200°C | -50°C - +200°C |
| Cooling Capacity @ 20°C | 1000 W | 1000 W |
| Work Area (D x W x L) mm | 200 x 208.5 x 104.2 | 200 x 208.5 x 104.2 |
| Bath Volume | 12 | 12 |
| Weight (kg) | 62 | 62 |

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.