



Product:

HAAKE Poly Lab OS
Docking Station

Manufactured By:

Thermo Fisher Scientific

Dimensions (L x W x H)	40 x 50 x 103 cm
Weight	50 Kg
Electrical Power	from RheoDrive base unit
Heating Zones	4 standard (more optional)
Electrical Heating/Air Cooling	3 sections
Max. Temperature	450°C *
Temperature Resolution	0.1°C
Max. Torque	300 Nm *
Max. Speed	200 min ⁻¹ *

* can be limited by the mixer / extruder itself

Supported Mixer & Extruder

Model	Year	Coupling
RheoMix 600 / 610 540	= 1990	557-1126
RheoMix 600 / 610 / 540	< 1990	557-1143
RheoMix 3000 / 3010	= 1990	557-1125
RheoMix 3000 / 3010	< 1990	557-1143
RheoMex 102 / 104 / 252 / 302	= 1990	557-1125
RheoMex 102 / 104 / 252 / 302	< 1990	557-1128
Conical Twin Screw Extruder TW100	< 1990	557-1141
Other Models	-	On request

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.