



Product:

HAAKE PolyLab Sheet and Ribbon Haul-Off

Manufactured By:

Thermo Fisher Scientific

Sheet and Ribbon Haul-off System

Roll diameter	110 mm
Roll width	200 mm
Roll surface	chrome – highly polished
Max. temperature	120°C or 200°C
Take off speed	0 - 10 m/min
Main supply	230V +/-10 % 50 - 60Hz
Max. current	6 A
Dimensions (L/W/H)	570/1130/720 mm
Weight	165 kg

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.