



Technical Specifications

Product:

HAAKE Sample Sizer

Manufactured By:

Thermo Fisher Scientific

Dimensions L x W x H	1200 x 960 x 1800 mm
Weight	450 Kg
Max. haul-off speed	1,0 m / min
Max. cutting rate	60 cuts / min
Max. profile dimensions W x H	Ø 30 mm, 40 x 15 mm, 50 mm ²
Electrical connection	400 V / 50 Hz / 3 ph + N + PE / 35 A
Vacuum pump	3 KW
Water connections	4
Vacuum connections	4
Air pressure connections	1/4" , 6 bar
Water inlet	3/4" , 4 bar
Water outlet	Ø 75 mm
Noise level	75 dB (A)

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.