

Product:

Full Notch Creep Test FNCT

Manufactured By:

Optical Control Systems

Load range	4 - 6 GPa on samples 10 x 10 x 100 mm (6 to 9 GPa on samples 6 x 6 x 90 mm optional)
Force	Resolution Indefinite, approved for 0,1 N, individual calibration by sample and position (customised) Calibration accuracy: better than +/- 1%
Fluid	Volume approx.: 55l Level control: stainless steel float sensors and solenoid valves Temperature range: RT to 95 °C Accuracy: 1 °C
Input pressure range for demin water supply	0,2 - 8 bar (3 - 116 psi)
Device interface	Ethernet 10/100 M Base T
Remote control	WEB interface
Size dimension	(l, w, h) 143 x 81 x 113 cm Weight approx. 510 kg
Power supply	230 V AC, 50/60 Hz
Temperature	10 - 40 °C

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.