

Product:

Solid Freeze / Thaw Chamber

Manufactured By:

Schleibinger Gerate

Outer Vessel size	220 mm
Inner Vessel size	Standard according to Standard SN 670 321a (2001)
Water coolers (No. of)	Two (2)
Refrigeration power	0.13 kW
Heating power	1.8 kW
Temperature Range	-10°C - +100°C
Temperature Deviation	< 0.02°C
LVDT Gauge	Better than 3 micron

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.