

Product:

Thin Layer Shrinkage Measurement
System

Manufactured By:

Schleibinger Gerate

Measurement range	2 * 4 mm
Specimen size	Max 30 x 25 cm
Resolution	0,1 μ m nm, class 1
Laser-spot diameter	0,8 mm
Security note	Laser power 1 mW at 625

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.