



Product:

Viskomat NT

Manufactured By:

Schleibinger Gerate

Stroke Range	5 mm
Sample Volume	370 ml
Speed	± 0.0001 to ± 250 rpm
Torque	0 to ± 250 Nmm
Control Unit (W x D x H)	190 x 450 x 350 mm
Viskomat (W x D x H)	300 x 670 x 50 mm
Weight	20 Kg

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.