



Product:

Shaking Water Bath

Manufactured By:

Thermo Fisher Scientific

Working Temperature Range: with tap water cooling	22°C - 90°C 20°C - 90°C
Temperature Accuracy	+/- 0.2°C
Heater Capacity (230V)	1.5 kW
Type of control	PID
Over temperature protection	Adjustable
Bath opening (W x L x D)	50 x 30 x 20 cm
Bath Volume	8 - 26 Litres
Shaking amplitude	15 mm
Shaking frequency (cont. adj.)	20 - 200 mm
Dimensions (overall)	65 x 34 x 26 cm
Net weight	18 Kg
Total Wattage (230V)	1600 VA

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.