

Product:

SAHARA Heating Bath Circulators

Manufactured By:

Thermo Fisher Scientific

Transparent Acrylic Models

Controller Head	Bath	S6T*	S12T*	S19T*
SC100		+13°C - +80°C	+13°C - +80°C	+13°C - +80°C
SC150		+13°C - +80°C	+13°C - +80°C	+13°C - +80°C
AC150		-	+13°C - +80°C	+13°C - +80°C
AC200		-	+13°C - +80°C	+13°C - +80°C
Maximum bath volume (litres)		6	12	19
Work area (D x W x L) mm		150 x 138 x 223	150 x 302 x 148.9	150 x 302 x 326.9
Weight (kg)		6.3	7.3	8.7

Accessories	Fluids, Test Tube Racks, Lifting Platform Tap Water Cooling Coil
--------------------	---

* Bath lids are NOT available for these models

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.